



Effective Date:
Monday, May 22, 2017

Test Updates

Immediate Action

In our continuing effort to provide you with the highest quality toxicology laboratory services available, we have compiled important changes regarding a number of tests we perform. Listed below are the types of changes that may be included in this notification, effective Monday, May 22, 2017

Test Changes - Tests that have had changes to the method/ CPT code, units of measurement, scope of analysis, reference comments, or specimen requirements.

Discontinued Tests - Tests being discontinued with alternate testing suggestions.

Please use this information to update your computer systems/records. These changes are important to ensure standardization of our mutual laboratory databases.

If you have any questions about the information contained in this notification, please call our Client Support Department at (866) 522-2206. Thank you for your continued support of NMS Labs and your assistance in implementing these changes.

The CPT Codes provided in this document are based on AMA guidelines and are for informational purposes only. NMS Labs does not assume responsibility for billing errors due to reliance on the CPT Codes listed in this document.



Effective Date:
Monday, May 22, 2017

Test Updates

Test Code	Test Name	Test Name	Method / CPT Code	Specimen Req.	Stability	Scope	Units	Reference Comments	Discontinue
5105U	Acetyl Fentanyl and Metabolite Confirmation, Urine (Forensic)	•			•	•			
52149U	Acetyl Fentanyl and Metabolite Confirmation, Urine (Forensic)	•			•	•			
0205U	Acetyl Fentanyl and Metabolite, Urine	•			•	•			



Effective Date:
Monday, May 22, 2017

Test Updates

Test Changes

5105U Acetyl Fentanyl and Metabolite Confirmation, Urine (Forensic)

Summary of Changes: Test Name was changed.
Stability was changed.
Scope of Analysis was changed.
Acetyl Norfentanyl was added.

Stability: Room Temperature: 30 day(s)
Refrigerated: 30 day(s)
Frozen (-20 °C): 30 day(s)

Scope of Analysis: LC-MS/MS (80362): Acetyl Fentanyl, Acetyl Norfentanyl
Method (CPT Code)

Compound Name	Units	Reference Comment
Acetyl Norfentanyl	ng/mL	

52149U Acetyl Fentanyl and Metabolite Confirmation, Urine (Forensic)

Summary of Changes: Test Name was changed.
Stability was changed.
Scope of Analysis was changed.
Acetyl Norfentanyl was added.

Stability: Room Temperature: 30 day(s)
Refrigerated: 30 day(s)
Frozen (-20 °C): 30 day(s)

Scope of Analysis: LC-MS/MS (80362): Acetyl Fentanyl, Acetyl Norfentanyl
Method (CPT Code)

Compound Name	Units	Reference Comment
Acetyl Norfentanyl	ng/mL	

0205U Acetyl Fentanyl and Metabolite, Urine

Summary of Changes: Test Name was changed.
Stability was changed.
Scope of Analysis was changed.
Acetyl Norfentanyl was added.

Stability: Room Temperature: 30 day(s)
Refrigerated: 30 day(s)
Frozen (-20 °C): 30 day(s)

Scope of Analysis: LC-MS/MS (80362): Acetyl Fentanyl, Acetyl Norfentanyl
Method (CPT Code)

Compound Name	Units	Reference Comment
Acetyl Norfentanyl	ng/mL	