



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Demo Report

Report Issued 05/14/2021 13:52
Last Report Issued 05/14/2021 13:52

88888
Clinical Example Report
Attn: Example Reports
200 Welsh Road
Horsham, PA 19044

Patient Name 1006TI-POS
Patient ID 1006TI-POS
Chain 20001570
Age Not Given **DOB** Not Given
Gender Not Given
Workorder 20001570
Received 07/22/2020 07:29

Sample ID 20001570-001
Matrix Tissue
Patient Name 1006TI-POS
Patient ID 1006TI-POS
Container Type Clear vial

Collect Dt/Tm Not Given
Source Not Given

Approx Vol/Weight Not Given

Receipt Notes None Entered

Analysis and Comments	Result	Units	Reporting Limit	Notes
-----------------------	--------	-------	-----------------	-------

1006TI Carbon Monoxide - Iron Ratio Profile, Tissue

Analysis by Gas Chromatography/Mass Spectrometry (GC/MS)

Carbon Monoxide	5.0	mcg/g	1.1	
-----------------	-----	-------	-----	--

Analysis by Inductively Coupled Plasma/Optical Emission Spectrometry (ICP/OES)

Iron	None Detected	mcg/g	20	
------	---------------	-------	----	--

Carboxyhemoglobin	None Detected	%Saturation		
-------------------	---------------	-------------	--	--

Synonym(s): COHb

This result is a calculated result based on the Iron and Carbon Monoxide analyses. The calculated Carboxyhemoglobin percent saturation underestimates the actual Carbon Monoxide percent saturation since the determination assumes that all Iron present is capable of binding Carbon Monoxide.

Results for sample 20001570-001 are continued on next page



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Sample ID 20001570-001
Matrix Tissue
Patient Name 1006TI-POS
Patient ID 1006TI-POS

Collect Dt/Tm Not Given
Source Not Given

Analysis and Comments	Result	Units	Reporting Limit	Notes
-----------------------	--------	-------	-----------------	-------

Specimens for elemental testing should be collected in certified metal-free containers. Elevated results for elemental testing may be caused by environmental contamination at the time of specimen collection and should be interpreted accordingly. It is recommended that unexpected elevated results be verified by testing another specimen.

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the US Food and Drug Administration.